

PICO Amplification Core Kit

The **PICO Amplification Core (AMC) Kit** (#PICO-000010) contents must be stored at 4°C or -20°C (see product labels for exact storage conditions). Find the date of expiry on the back side of the pouches.

Further information:

- The PICO AMC User Manual, the PICO Calculator, and the Safety Data Sheets are available at: www.actome.de/resources/downloads
- The AMULATOR can be accessed here: <https://amulator.actome.de/>

Notes before starting:

- The PICO AMC workflow is compatible with both 24- and 96-well QIAcuity Nanoplates.
- Perform all centrifugation steps with 400 rcf for 5 min at RT unless stated otherwise.
- **Important:** For absolute quantification (AQ), we recommend an antibody concentration of 5×10^{-10} M in the binding reaction. If AQ is *not* required, we recommend an antibody concentration of 4×10^{-11} M.
- Mixing by pipetting should be done slowly (30 times) to avoid generation of air bubbles.
- We recommend using at least four technical replicates for each sample.
- For the analysis of cultivated cells, we recommend using a total of 1 million cells in a concentration of 1×10^4 cells/ μ l as an input.
- We recommend using 3-fold dilution steps for high abundance proteins and 2-fold dilution steps for low abundance proteins.
- Sample preparation is dependent on the sample type used. Here, cultivated cells are described.
- During the analysis of protein interactions using a cross-linker is optional.

DAY 1

Preparation of biological material:

1. Prepare the following chemicals and buffers:

Additive C (5x stock)	BSA (5x stock)	EDTA-free Protease Inhibitor Cocktail (PIC), (25x stock) - not provided	PIC-PBS
add 500 μ l PBS	add 400 μ l PBS	dissolve 1 tablet of cOmplete Protease™ Inhibitor Cocktail (PIC) in 2 ml PBS	dissolve 1 PIC tablet in 500 μ l PBS
Lysis Buffer Stock (LB-Stock), (2x stock)	Lysis Buffer (LB)	Control Buffer (CB)	BS3 Stock Solution (BS3S), (20x stock) - optional
200 μ l Additive T 400 μ l Additive C 80 μ l PIC 200 μ l Additive L 120 μ l PBS	300 μ l LB-Stock 300 μ l PBS	250 μ l LB-Stock 100 μ l BSA 150 μ l PBS	dissolve 2 mg BS3 in 37.8 μ l DMSO
			BS3 Working Solution (BS3W) - optional
			6 μ l BS3S 114 μ l PBS

2. Harvest the cells with the preferred method and wash them in **~5 ml PIC-PBS**. Centrifuge and discard supernatant.
3. Resuspend the cells in **1 ml PIC-PBS** and transfer them to a 1.5 ml reaction tube. Centrifuge and discard supernatant.
4. Repeat the washing step with **1 ml PIC-PBS**, centrifuge, and discard supernatant. If cross-linking is required use **PBS** instead and continue with step 5, otherwise skip ahead to step 8.

5. Resuspend the cells in **1 ml PBS**, count and transfer 2×10^6 cells into a new 1.5 ml reaction tube. Centrifuge and discard supernatant.
6. Add **100 μ l BS3W**, mix by pipetting, and incubate for 30 min at RT.
7. Add **900 μ l PIC-PBS**, centrifuge, and discard supernatant.
8. Resuspend the cells in **1 ml PIC-PBS**, count, and transfer 1×10^6 cells into a new 1.5 ml reaction tube. Centrifuge and discard supernatant carefully.
9. Resuspend the cells in **100 μ l LB** and vortex for 10 s. Lyse the cells for 3 h at 4°C.
10. Sonicate the lysate for 5 min at full power in an ultrasonic bath at RT.
11. Transfer the lysate into a QIAshredder spin column. Centrifuge ($\sim 20,000$ rcf, 2 min) and transfer the flow-through to a new 1.5 ml reaction tube.

Binding reaction:

12. Calculate the volume of antibody stocks, LB-Stock, and PBS required for the antibody mix (ABX) using the **PICO Calculator** (step 2 in the 'PICO Calculator' tab). Make sure the correct Nanoplate is selected in cell E7.
13. Set up the binding reaction in a 96-well PCR microplate (see example layout for a 24-well PICO assay in **Figure 1**). For the **NTC** add **4 μ l LB** to the A1 well. For the **ABC** combine **2 μ l CB** with **2 μ l ABX** (prepare three technical replicates). For the **samples** combine **2 μ l biological sample** with **2 μ l ABX** (prepare four technical replicates).
14. Seal the plate with adhesive foil and sonicate for 1 min at full power in an ultrasonic bath at RT. Centrifuge the plate ($\sim 1,000$ rcf, 30 s) and incubate at 4°C overnight.

Pre-dilution and digital PCR:

15. Prepare the **Master Mix** for 24 or 96 samples, vortex for 10 s, and quick spin ($\sim 1,000$ rcf, 5 s).

Important: *If less than four PICO Probes are used, replace the volume of not used PICO Probe(s) with ultrapure water.*

Master Mix		
Reagents	24-well plate	96-well plate
Ultrapure water	606 μ l	834 μ l
QIAcuity Probe Master Mix	284 μ l	390 μ l
PICO Probe (P8, BL, N6, or O7)	45 μ l	62 μ l
PICO Probe (P8, BL, N6, or O7)	45 μ l	62 μ l
PICO Probe (P8, BL, N6, or O7)	45 μ l	62 μ l
PICO Probe (P8, BL, N6, or O7)	45 μ l	62 μ l
Coupling dPCR Mix	36 μ l	62 μ l

16. Calculate the required dilutions to target a **lambda of 0.15** in step 3 of the **PICO Calculator**.
17. For a 24-well PICO assay prepare a new 96-well plate for the dilution steps. Add the calculated amount of **PBS** to columns 1 - 3 and 5 - 7, and **41 μ l Master Mix** to columns 10 - 12 (**Figure 2**). For a 96-well PICO assay prepare three additional 96-well plates. Add the calculated amount of **PBS** to Dilution plate #1 and #2 and add 13 μ l Master Mix to the third 96-well plate (**Figure 3**).
18. Remove the adhesive foil of the plate containing the binding reaction and add **36 μ l PBS** to all wells (reflects the '10x pre-dilution of binding reaction' in **Figure 2 and 3**). Mix thoroughly by pipetting.
19. Transfer the predetermined volume (usually 1 μ l) from each pre-diluted sample to the corresponding wells of the dilution plate and mix by pipetting (DS1/Dilution plate #1). Repeat the dilution once more (DS2/Dilution plate #2) and finally transfer **1 μ l** to the wells/plate containing the **Master Mix (Figure 2 and 3)**. Mix by pipetting.

20. Transfer **40 µl or 12 µl Master Mix**, containing the diluted samples, to a 24- or 96-well QIAcuity Nanoplate respectively. Seal the plate according to the QIAcuity user manual and run a dPCR using the following parameters:

Priming			
QIAGEN Standard Priming Profile			
PCR conditions			
Step	Temperature	Time	
Hot-start	95°C	2 min	Cycle 40 times
Denaturing	95°C	15 s	
Annealing	58°C	30 s	

Imaging conditions			
PICO Probe	QIAcuity channel	Integration time	Gain
P8 Probe	FAM, green channel	500 ms	6
BL Probe	HEX, yellow channel	400 ms	6
N6 Probe	TAMRA, orange channel	400 ms	6
O7 Probe	ROX, red channel	300 ms	4

21. Instructions for the data analysis with AMULATOR can be found in the user manual.

Binding reaction




Figure 1: Plate layout of an example binding reaction for a 24-well PICO assay. In this example five samples, with four technical replicates each, were used. For the ABC three technical replicates are sufficient. ABC - antibody control; NTC - non template control




Figure 2: Exemplary setup to perform the dPCR pre-dilution for a 24-well PICO assay. The ‘10x pre-dilution of binding reaction’ is prepared in the plate containing the binding reaction.




Figure 3: Exemplary setup to perform the dPCR pre-dilution for a 96-well PICO assay. The '10x pre-dilution of binding reaction' is prepared in the plate containing the binding reaction.

Scan the QR code for the user manual:

